



act your doctor or physical
ree of discomfort is normal. Talk to
6 to 12 weeks, unless otherwise hysical therapist can instruct you is on and health of your elbows.
ities, as well as sports and other
n will help with the healing rance and improve resistance to
e painful conditions caused by ng the same motions again and
med under your doctor's help you meet your rehabilitation
ndylitis

AAOS does not endorse any treatments, procedures, products, or physicians referenced herein. This information is provided as an educational service and is not intended to serve as medical advice. Anyone seeking specific orthopaedic advice or assistance should consult his or her orthopaedic surgeon.





1. Wrist Extension Stretch

Repetitions

Equipment needed: None

5 reps, 4x a day

Additional instructions: This stretch should be done throughout the day, especially before activity. After recovery, this stretch should be included as part of a warm-up to activities that involve gripping, such as gardening, tennis, and golf.


Days per week

5 to 7

Step-by-step directions

- Straighten your arm and bend your wrist back as if signaling someone to "stop."
- Use your opposite hand to apply gentle pressure across the palm and pull it toward you until you feel a stretch on the inside of your forearm.
- Hold the stretch for 15 seconds.
- Repeat 5 times, then perform this stretch on the other arm.

Tip Do not lock your elbow.



2. Wrist Flexion Stretch

Repetitions

Equipment needed: None

5 reps, 4x a day **Additional instructions:** This stretch should be done throughout the day, especially before activity. After recovery, this stretch should be included as part of a warm-up to activities that involve gripping, such as gardening, tennis, and golf.


Days per week

5 to 7

Step-by-step directions

- Straighten your arm with your palm facing down and bend your wrist so that your fingers point down.
- Gently pull your hand toward your body until you feel a stretch on the outside of your forearm.
- Hold the stretch for 15 seconds.
- Repeat 5 times, then perform this stretch on the other arm.

Tip Do not lock your elbow.





3. Wrist Extension (Strengthening)

Repetitions

Equipment needed: Dumbbell hand weights (1 lb., 2 lbs., 3 lbs.)

30 reps, 1x a day

Additional instructions: This exercise should be performed in stages. Begin each stage with no weight. When you are able to complete 30 repetitions on 2 consecutive days with no increase in pain, move forward in the program by increasing weight (begin with 1 lb., advance to 2 lbs., end with 3 lbs.).

Days per week 5 to 7


Stage 1: Bend your elbow to 90 degrees and support your forearm on a table with your wrist placed at the edge. Early in rehabilitation for painful lateral epicondylitis, just doing the lowering down portion of this exercise using a light (1 to 3 lb.) dumbbell can help resolve acute pain before progressing to doing the complete motion once there is less pain. Use the opposite hand to help push the wrist up as far as possible. Then, without assistance, slowly lower the weight until the wrist is flexed again.

Stage 2: Straighten your elbow slightly. Continue to support your arm on the table.

Stage 3: Fully straighten your elbow and lift your arm so that it is no longer supported by the table.

Step-by-step directions to be followed for each stage


- With palm down, bend your wrist up as far as possible.
- Hold up for 1 count, then slowly lower 3 counts.
- Begin without using a weight and increase the repetitions until you can complete 30.
- When you can perform 30 repetitions on 2 consecutive days without increasing pain, begin performing the exercise using a 1 lb. weight.
- Follow the same steps above to continue to increase repetitions and weight until you are using a 3 lb. weight and can perform 30 repetitions on 2 consecutive days without increasing pain.



Stage 3 using 1 lb. weight

Stage 1 using 1 lb. weight

Tip Do not let the weight pull your hand down too quickly.







4. Wrist Flexion (Strengthening)

Repetitions

Equipment needed: Dumbbell hand weights (1 lb., 2 lbs., 3 lbs.)

30 reps, 1x a day

Additional instructions: This exercise should be performed in stages. Begin each stage with no weight. When you are able to complete 30 repetitions on 2 consecutive days with no increase in pain, move forward in the program by increasing weight (begin with 1 lb., advance to 2 lbs., end with 3 lbs.).

Days per week

5 to 7

Stage 1: Bend your elbow to 90 degrees and support your forearm on a table with your wrist placed at the edge. Early in rehabilitation for painful medial epicondylitis, just doing the lowering down portion of this exercise using a light (1 to 3 lb.) dumbbell can help resolve acute pain before progressing to doing the complete motion once there is less pain. Use the opposite hand to help push the wrist up as far as possible. Then, without assistance, slowly lower the weight until the wrist is flexed again.


Stage 2: Straighten your elbow slightly. Continue to support your arm on the table.

Stage 3: Fully straighten your elbow and lift your arm so that it is no longer supported by the table.


Step-by-step directions to be followed for each stage

- With palm up, bend your wrist up as far as possible.
- Hold up for 1 count, then slowly lower 3 counts.
- Begin without using a weight and increase the repetitions until you can complete 30.
- When you can perform 30 repetitions on 2 consecutive days without increasing pain, begin performing the exercise using a 1 lb. weight.
- Follow the same steps above to continue to increase repetitions and weight until you are using a 3 lb. weight and can perform 30 repetitions on 2 consecutive days without increasing pain.





Stage 1 using 1 lb. weight



Stage 3 using 1 lb. weight





5. Forearm Supination & Pronation (Strengthening)

Repetitions

Equipment needed: Dumbbell hand weights (1 lb., 2 lbs., 3 lbs.)

30 reps, 1x a day

Additional instructions: This exercise should be performed in stages. Begin each stage with no weight. When you are able to complete 30 repetitions on 2 consecutive days with no increase in pain, move forward in the program by increasing weight (begin with 1 lb., advance to 2 lbs., end with 3 lbs.).

Days per week

5 to 7


Stage 1: Bend your elbow to 90 degrees and support your forearm on a table with your wrist placed at the edge.

Stage 2: Straighten your elbow slightly. Continue to support your arm on the table.


Stage 3: Fully straighten your elbow and lift your arm so that it is no longer supported by the table.

Step-by-step directions to be followed for each stage

- Begin with palm facing the side. Slowly turn the palm facing up.
- Slowly return to the start position, then slowly turn the palm down.
- Slowly return to start position. This completes one repetition.
- Begin without using a weight and increase the repetitions until you can complete 30.
- When you can perform 30 repetitions on 2 consecutive days without increasing pain, begin performing the exercise using a 1 lb. weight.
- Follow the same steps above to continue to increase repetitions and weight until you are using a 3 lb. weight and can perform 30 repetitions on 2 consecutive days without increasing pain.



Stage 1 using 1 lb. weight



Stage 3 using 1 lb. weight

Tip

When using a dumbbell, try to let the weight pull your forearm in either direction as far as possible.






7. Stress Ball Squeeze_

Repetitions

Equipment needed: Rubber stress ball

10 reps, 1x day

Additional instructions: This exercise should be performed after completing the above staged strengthening exercises. Your arm and elbow position for this exercise should match the stage you are completing.





8. Finger Stretch

Repetitions

Equipment needed: Elastic band

10 reps, 1x day

Additional instructions: This exercise should be performed after completing the above staged strengthening exercises. Your arm and elbow position for this exercise should match the stage you are completing.

Days per week

5 to 7



